# milenia biotec

## HybriDetect 2T

**Universal lateral flow dipstick for simultaneous detection of two different analytes (proteins, genomic amplicates) labelled with FITC, biotin and digoxigenin; development platform**  
English: Page 1-10

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E-Mail: [info@milenia-biotec.de](mailto:info@milenia-biotec.de)

Web: [www.milenia-biotec.de](http://www.milenia-biotec.de)


**Note:** Significant changes are indicated by dotted lines in the margin. In the end of the IFU you will find a table with causes of changes.

## Table of Contents

Explanation of Symbols .....	3
Warnings and Precautions .....	3
Intended Use .....	3
Materials Supplied, Storage and Stability .....	4
Materials Required (not included) .....	4
Necessary Development Work - Development Platform .....	4
Method .....	5
Testprinciple – Genamplicon Detection .....	6
Testprinciple – Protein Detection .....	7
Control Band Dipstick .....	7
Assay Performance “PCR Products” .....	8
Assay Performance “RPA Products” .....	8
Interpretation of Results .....	8
Trouble Shooting “PCR” .....	9
Assay Sensitivity “PCR” .....	9
Literature-References .....	9
Additional available Products .....	10

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## Warnings and Precautions

- Store reagents at **2°C to 8°C** in original packaging
- Before use, bring all reagents to room temperature (18 - 28° C).
- Do not use reagents after the expiry date has reached!
- Protect dipsticks from humidity, Container must always be closed
- Touch and label only the foil-covered areas of the dipsticks (labeling area).
- The disposal of waste materials must be carried out according to current local regulations.
- The assay buffer contains an anti-microbial reagent; therefore avoid contact with skin and/or mucous membranes.
- For professional users

## Intended Use

HybriDetect 2T is a universal lateral flow dipstick for simultaneous detection of two different analytes (proteins, genomic amplicates) labelled with FITC, biotin and digoxigenin. It is a development platform. The test is for research use only, not for diagnostic purpose.

## Materials Supplied, Storage and Stability

Component	Cat.-No.	Content	Preparation	Storage	Shelf Life
<b>HybriDetect 2T Dipsticks:</b> Membrane coated with biotin-ligand and polyclonal (goat) digoxigenin antibody; Conjugate: polyclonal (goat) anti-FITC antibody labelled gold particles	MGDS2A	2 x 50 tests	Ready-to use	2 - 8°C; Containers must always be closed (moisture protection)	until expiry date
<b>HybriDetect 2T Assay Buffer:</b> Citrate-Phosphate Buffer	MGCBB	2 vials à 10 mL	Ready-to use	2 - 8°C	until expiry date

Material safety data sheets are available on request or from the company website ([www.milenia-biotec.de](http://www.milenia-biotec.de))

## Materials Required (not included)

- Pipets
- Pipet tips (containing protective filters for PCR)
- Reaction tubes or 96-well microtiter plate

## Necessary Development Work - Development Platform

Development of two solutions (A and B). Each solution must contain two different labeled detectors for the respective analyte.

The following conditions are necessary:

1) **Solution A:** Detectors must be labeled with:

- FITC (fluorescein isothiocyanate) or FAM (Carboxyfluorescein)
- Biotin

2) **Solution B:** Detectors must be labeled with:

- FITC (fluorescein isothiocyanate) or FAM (Carboxyfluorescein)
- Digoxigenin

3) Use about 100 µL fluid (sample material and analyte-specific solutions, A and B) for the assay procedure.

4) The provided buffers (MGCBB) may be used as a basic buffers for these analyte-specific solutions.

Application example:

20 µL sample material and 80 µL analyte-specific solutions with 5 min incubation time.

### **Notice:**

Volumes, analyte-specific solutions, and incubation time are part of the individual test development.

A basic procedure for the detection of genomic amplicons is explained on page 7: "Assay performance PCR products." One amplification product should be biotinylated and the specific hybridization probe should be FITC labeled. The second amplification product should be labelled with digoxigenin and the specific hybridization probe should be FITC labeled. Various amplification procedures (**PCR** or isothermal amplifications like **LAMP** or **RPA**) can be used.


## Method

Milenia HybriDetect 2T is a ready-to-use, universal test strip (dipstick), which is based on the lateral flow technology using gold particles. The dipstick is designed to develop qualitative or quantitative rapid test systems for simultaneous detection of two different analytes such as proteins, antibodies, or gene amplifications. The user needs to develop two analyte-specific solutions. Solution A: Contains a first detector (e.g. antibody, antigen, specific probe) labeled with FITC and a second one (e.g. antibody, primer) labeled with biotin. Solution B: Contains a first detector (e.g. antibody, antigen, specific probe) labeled with FITC, too and a second one (e.g. antibody, primer) labeled with digoxigenin (see “Common Testprinciple” on page 5).


The sample to be determined is mixed with the developed analyte-specific solutions and the dipstick is placed into this solution.

The analyte A-complex labeled with FITC and biotin binds first to the gold-labeled FITC-specific antibodies in the sample application area of the dipstick. The analyte B-complex (labeled with FITC and digoxigenin) also binds to this area of the dipstick. Capillary forces cause the gold complexes A and B to diffuse across the analytical membrane. Only the analyte captured gold particles will bind when they overflow the immobilized biotin-ligand molecules at the respective test band (test band A- analyt A, test band B- analyt B) and generate a red-blue band over the time. Not-captured gold particles flow over the control band and will be fixed there by species-specific antibodies. With increasing incubation time, the formation of an intensely colored control band appears.

## Testprinciple – Genamplicon Detection

## Testprinciple – Protein Detection


## Control Band Dipstick

**In any case, the control band must be visible!**

It is a control function and cannot be used to assess the test band. If the control band is not visible after the incubation period, the result is invalid! The test must be repeated with a new dipstick!

## Assay Performance “PCR Products”

1. Take the required number of dipsticks out of the container and mark them.
2. For each sample to be analyzed pipet **100 µL** HybriDetect Assay Buffer or individual developed buffer into a reaction tube or a well of a microtiter plate.
3. Pipet **5 - 10 µl** of the hybridization product directly on the sample application area, or **alternatively**: add **5 - 10 µl** of the hybridization product into the solution of the reaction tube / well.
4. Place the dipsticks with the sample application area into the solution and incubate them e.g. for **5 - 15 minutes** in an upright position.
5. In the end of incubation period, remove dipsticks from assay solution and interpret test results immediately.

### Note:

- If a higher analytical sensitivity is required, it could be helpful to increase the volume of the PCR product.
- Volumes, analyte-specific solutions, and incubation time are always part of the individual test development.

## Assay Performance “RPA Products”


1. Take the required number of dipsticks out of the container and mark them.
2. For each sample to be analyzed pipet **80 µL** HybriDetect Assay Buffer (MGCBB) into a reaction tube or a well of a microtiter plate.
3. After RPA reaction dilute the RPA reaction mix 1/50 or 1/25 (for a higher sensitivity) with HybriDetect Assay Buffer (MGCBB) and load 10 µl directly on the sample application area.
4. Place the loaded dipsticks with the sample application area into the prepared solution (see 2.) and incubate them e.g. for 5 - 15 minutes in an upright position.

In the end of incubation period, remove dipsticks from assay solution and interpret test results

## Interpretation of Results

### Applying an internal PCR control:

While the biotin-labeled **amplicon A** is mainly used for the detection of a specific target sequence the digoxigenin-labeled **amplicon B** is provided as an internal PCR amplification control.

Test band A	Test band B	Control band	Interpretation
positive	positive	positive	<ul style="list-style-type: none"> <li>control band is clearly visible, test run is valid</li> <li><b>amplicon A is detected (positive)</b></li> <li>internal PCR control is positive (amplicon B)</li> </ul>
negative	positive	positive	<ul style="list-style-type: none"> <li>control band is clearly visible, test run is valid</li> <li><b>amplicon A is not detected (negative)</b></li> <li>internal PCR control is positive (amplicon B)</li> </ul>
positive	negative	negative	<ul style="list-style-type: none"> <li>test runs were not valid because control band did not appear</li> </ul>
negative	positive		
negative	negative		
negative	negative	positive	<ul style="list-style-type: none"> <li>rapid test performed properly because control band is clearly visible</li> <li>PCR amplification did not work due (amplicon B is not detectable )</li> </ul>
positive			

## Trouble Shooting “PCR”

Problem	Possible cause(s)	Recommendation
Control band is not visible.	<ul style="list-style-type: none"> <li>a) wrong or destroyed assay buffer</li> <li>b) expiration date of dipsticks is exceeded</li> <li>c) wrong storage conditions of dipsticks</li> </ul>	apply new (fresh) chemicals
Negative result with dipstick but clearly visible band in agarose gel	<ul style="list-style-type: none"> <li>a) detection of an unspecific PCR product in agarose gel</li> <li>b) hybridization was not successful</li> </ul>	Check identity of PCR product by Southern blotting or sequence analysis Check conditions of hybridization.
Mineral oil	<ul style="list-style-type: none"> <li>a) mineral oil affects flow characteristics of the assay</li> <li>b) development of band might be hampered</li> </ul>	Remove PCR product very slowly from the bottom of the reaction vial.

## Assay Sensitivity “PCR”

Analytical sensitivity of the Milenia HybriDetect 2T is equivalent to agarose gel electrophoresis and subsequent staining with ethidiumbromide. Independent on the size of the PCR product and the number of amplification cycles as low as 5 pg DNA could be detected.

## Literature-References

Dai T, et al, *Frontiers in Microbiology* (2019); 10:1884.  
 Yang Y, et al, *Virology Journal* (2017), 14(1), 1–10.  
 Qi Y, et al, *PLoS ONE* (2018), 13(11).  
 Ahmed F A, et al, *Scientific Reports* (2018), 8(1).  
 Wu L, *Parasites & Vectors* (2019). 12(1),

## Additional available Products

Product Name	Order No.	Content	Description
Milenia HybriDetect	MGHD 1	100 tests	Dipsticks with <b>one</b> test band

For any questions and suggestions, please contact us:

**Telefon:** +49 641 – 94 88 83 – 0

**Fax:** +49 641 – 94 88 83 – 80

**Web:** [www.milenia-biotec.de](http://www.milenia-biotec.de)

milenia biotec

## HybriDetect 2T

Universal lateral flow dipstick for simultaneous detection of two different analytes (proteins, genomic amplicates) labelled with FITC, biotin and digoxigenin; development platform  
English: Page 1-10

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Milenia Biotec GmbH · Versailler Str. 1 · 35394 Gießen, Germany