Datasheet

Cat# NBGN-100132 Version# RN5.9



Product Name	NebuChem™ Exosome Isolation Magnetic Beads
Size	1ml; Bulk
Description	NebuChem™ Exosome Isolation Magnetic Beads (CAT#NBGN-100132) is used for
	one-step exosome isolation from serum, plasma, saliva, or other bodily fluids.
Components	High-Efficient Exosome Binding Beads, 1ml (CAT#NBGN-100132A)
	5X Binding Buffer, 2.5ml (CAT#NBGN-100132B)
	10X Washing Buffer, 7.5ml (CAT#NBGN-100132C)
	Elution Buffer, 2ml (CAT#NBGN-100132D)
Principle	STEP1: CentrifugationHost cells are separated and removed.
	STEP2: Bead CaptureMagnetic beads are utilized to capture and trap exosomes.
	STEP3: WashingThe beads undergo washing to eliminate nonspecific binding.
Protocol	Materials Needed But Not Provided
	1) Sterile 1.5ml Centrifuge Tubes
	2) Centrifuge
	3) Micropipettes
	4) RNase&DNase-Free Water
	5) Magnetic Rack
	Step A: Reagents Preparation
	1) Before pipetting, vortex to completely resuspend the magnetic beads.
	2) Do not allow the magnetic beads to sit for more than 2 minutes before
	dispensing.
	3) Dilute all the stock buffer with ultra-pure water to 1X working buffer
	Step B: Sample Preparation
	NebuChem™ Exosome Isolation Magnetic Beads is optimized for processing
	$\sim\!250\mu L$ of serum, plasma, saliva, or other bodily fluids.
	Step C: Isolation & Purification
	1) Transfer 250μL of biological samples (serum, or plasma) to a 1.5 ml centrifuge
	tube.
	2) Centrifuge at 12,000 x g for 10 minutes and carefully transfer the supernatant
	to a new 1.5 ml centrifuge tube without disturbing the pellet.
	3) Add 250 μL of 1x Binding Buffer and 20μL of Magnetic Beads.
	Note: Vigorously shake the bottle until the magnetic beads become

Datasheet

Cat# NBGN-100132 Version# RN5.9



homogeneous before dispensing. Do not allow the beads to sit for more than 2 minutes before dispensing. Resuspend the magnetic beads every 2 minutes.

- 4) Mix the beads by gently pipetting up and down to mix for 20-25 times.
- 5) Separate the beads with a magnet and remove the supernatant.
- 6) Add 500 μ l of Washing Buffer and wash the beads by pipetting up and down 5 -10 times.
- 7) Separate the beads with a magnet and remove the supernatant.
- 8) Repeat steps 6-7 to wash the beads three times.
- 9) Suspend the beads with 50 μ L of Elution Buffer and pipet up and down for 20 times.
- 10) Separate the beads with a magnet and transfer the supernatant containing exosome lysates to a new centrifuge tube.to a fresh tube.
- 11) The exosome is ready for downstream applications.

Storage

2 - 8°C; Freezing is prohibited